



STERIS FREEZE DRYER
cGMP COMPLIANT
SIP, CIP

NEW

Interior



Fourteen and One Half (14.5)
Square Feet of Useable Shelf
Area, Stoppering System
For Vials, 20 Kg Ice Condenser

VIAL CAPACITY (Approx)

5 ml	2,700
10 ml	2,160
20 ml	1,330

SPECIFICATIONS

Number of Useable Shelves:	5
Dimensions of Shelves (mm):	450 wide, 600 deep
Temperature Range, Shelves:	- 55 C to 70 C
Final Condenser Temperature:	- 75 C
Condenser Capacity:	20 Kg
Thermal Medium:	Silicone Oil
Chamber Design:	Rectangular
Pressure Range:	Amb to 30 psig
Vacuum Range	Amb to 10 ⁻⁴ Bar
Instrumentation:	PC Interface to Allen Bradley PLC
Construction:	316 L Polished Stainless Steel, (Ra less than 1 m/m)
Electrical:	460/3/60
Connected Load:	36 kVA
Maximum Amperes:	50

American Lyophilizer, Inc
668 Stony Hill Road
Yardley, PA 19067
e-mail: lyophi@comcast.net
Telephone: 610 999 4151